
# **Praktikumsskript**

## **Masterpraktikum**

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**Wintersemester 2025/26**

# Inhalt

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

## 5. Synthesis of Precursor: 1,2-HOPO(Bn)-thiazolidine

### Monday - First Day:

1,2-HOPOBn (\_\_\_\_\_ mg, 0.50 mmol, 1.00 eq) was dried overnight *in vacuo*.

### Tuesday - Second Day:

The dried 1,2-HOPOBn was suspended in dry toluene (2 mL). Oxalyl chloride (\_\_\_\_\_  $\mu\text{L}$ , 0.90 mmol, 1.8 eq) was added at room temperature under an  $\text{N}_2$  atmosphere. Under stirring, 2-3 drops of DMF were added every 1 hour until no further gas evolution was observed. The solvent and excess oxalyl chloride were carefully removed *in vacuo* overnight. The acid chloride was then used without further purification immediately in the next step.

### Wednesday - Third Day:

The 1,2-HOPOBn acid chloride was dissolved in dichloromethane (5 mL) and added dropwise to a vigorously stirred solution of mercaptothiazolidine (\_\_\_\_\_ mg, 0.60 mmol, 1.20 eq) and potassium carbonate (\_\_\_\_\_ mg, 2.50 mmol, 5.00 eq) in water (1.5 mL) while cooling in an ice bath. After complete addition, the bright yellow mixture was warmed to room temperature and stirred overnight.

### Thursday - Fourth Day:

The layers were separated and the aqueous layer extracted with dichloromethane ( $3 \times 10$  mL). The combined organic phases were dried over sodium sulfate and filtered. The crude product was analyzed by TLC {DCM & MeOH} and then evaporated under reduced pressure using a rotary evaporator.


The samples, a yellow oil collected from all students, will be combined and purified by the assistant using column chromatography.

## 6. CAS (chrome azurol S)-based overlay assay

**Background:** Siderophores are chelators secreted by prokaryotes to mobilize poorly soluble iron and facilitate its subsequent uptake. Related chelators are known for other metals, e.g. chalkophores (copper) or zincophores (zinc).

In this part of the practical course, we will screen three different bacterial strains: *Escherichia coli* JM109, *Pseudomonas allopuntida* KT2440, and *Methylobacterium extorquens* PA1 for siderophore production.

We will take advantage of the metal-sensitive dye chrome azurol S (CAS) for the screening. We will grow the three strains on a medium lacking iron, thus stimulating siderophore production and release, and visualize siderophore release using CAS (Figure 1). CAS is a metal sensitive-dye and appears blue if iron is bound. If iron is released, the color changes to orange-yellow.

**Figure 1:** The overlay CAS assay in a nutshell. Cultures grown on a solid medium are overlaid with CAS-containing soft agar. Siderophores in the extracellular space mobilize iron from the CAS and cause a color change from blue to orange.

### Monday - First Day:

You will find all the needed materials on the bench assigned to your group. Details regarding the composition of the used media reagents can be found in the references listed at the end of this document.

We will start with a little background session about the essential elements of life, and how iron fits in here. Next, we will do a microbiology 101 crash course concerning cultivation and sterile work.

### Procedure:

1. Inspect the three plates we provided you with (one for each strain); how do the strains differ in terms of morphology (color, colony shape, etc.)?

Note: We will use this information on Thursday to assess growth and check for potential contamination.

2. Clean the surface of your bench with Bacillol® (a mixture of 1-, 2-propanol and ethanol).

Note: Don't spray the product on the bench (in microbiology we always want to avoid aerosol formation); spray it on a paper towel and wipe the surface.

3. Turn on the Bunsen burner and open the air hole fully. The roaring blue flame creates a zone of protection that allows you to work under sterile conditions. The rising air pushes microbes and particles up and away from your workspace. The resulting sterile workspace has roughly a radius of ~ 50 cm (Figure 2).

**Figure 2:** The roaring blue flame of the Bunsen burner creates a sterile workspace suitable for microbiological work.

Picture credit: neosynbio.com

4. Open one of the plates, pick up the biomass of a single colony with an inoculation loop, and streak it out on a fresh agar plate according to the four-streak technique.

Note: The idea is to distribute and “dilute” biomass. The latter leads ideally to single colonies, which facilitates checking for contaminants.

Prepared plates are incubated at 28 °C (*M. extorquens* PA1 and *P. alloputida* KT2440) and 37°C (*E. coli* JM109) until Thursday.

### **Wednesday - Third Day:**

#### Procedure:

1. Preheated (50 °C) CAS-containing soft agar is poured over the incubated plates and the plate surface completely covered. The soft agar layer should have a thickness of approximately 0.5 cm.
2. We inspect the plates for potential color changes after ~1 hour.

### **Thursday - Fourth Day:**

#### Procedure:


1. We re-visit our plates. After one day the colour change becomes more pronounced and plates can be interpreted more easily.

Note: Do you see differences regarding the intensity/shade? Check Pérez-Miranda et al., 2007. What might these tell us about the present siderophores?

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## Literatur

- Louden BC, Haarmann D, Lynne AM. 2011. Use of blue agar CAS assay for siderophore detection. *J Microbiol Biol Educ* 12:51–53.
- Schwyn B, Neilands JB. 1987. Universal chemical assay for the detection and determination of siderophores. *Anal Biochem* 160:47–56.
- Xu J, Churchill DG, Botta M, Raymond KN. 2004. Gadolinium(III) 1,2-hydroxypyridonate-based complexes: toward MRI contrast agents of high relaxivity. *Inorg Chem* 43:5492–5494.
- Fan D, Fang Q. 2021. Siderophores for medical applications: Imaging, sensors, and therapeutics. *Int J Pharm* 597:120306.
- Pérez-Miranda S, Cabirol N, George-Téllez R, Zamudio-Rivera LS, Fernández FJ. 2007. O-CAS, a fast and universal method for siderophore detection. *J Microbiol Methods* 70:127–131.
- Delaney NF, Kaczmarek ME, Ward LM, Swanson PK, Lee M-C, Marx CJ. 2013. Development of an optimized medium, strain and high-throughput culturing methods for *Methylobacterium extorquens*. *PLoS One* 8:e62957.
- <https://schlenklinesurvivalguide.com/> (abgerufen am 09.10.2025)

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