MOPOS-Chemiker*innen

Praktikumsskript

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
2.7 Reversible Addition-Fragmentation Chain-Transfer (RAFT) Polymerization of *N*-acryloylmorpholine (NAM)

2.7.1 Introduction and Aim

Reversible-deactivation radical polymerization (RDRP) minimizes irreversible termination and transfer reactions during free-radical polymerization (FRP) by providing a dynamic equilibrium between propagating radicals and dormant species. RDRP enables uniform chain growth, resulting in polymers with controlled molecular weights and low dispersity, provided by a fast initiation compared to the rate of propagation. A major advantage of RDRP is the retained chain-end functionality which allows further chain extension, enabling the synthesis of complex architectures such as block, star, brush, and hyperbranched polymers. Some important examples of RDRP include nitroxide-mediated polymerization (NMP), atom-transfer radical polymerization (ATRP), and reversible addition–fragmentation chain-transfer (RAFT) polymerization.

RAFT is particularly versatile among these due to its compatibility with a wide range of monomers, to which the structure of the so-called RAFT agent needs to be adapted. More-activated monomers (MAMs) show greater reactivity with radicals and form more stable macro-radicals, thus requiring more active RAFT agents (dithioesters and trithiocarbonates). In contrast, the polymerization of less-activated monomers (LAMs) is better controlled using less active RAFT agents with slower fragmentation rates (xanthates and dithiocarbamates). The success of RAFT polymerization depends on the ratios of monomer to RAFT agent and of RAFT agent to initiator.

The aim of this experiment is the synthesis of hydrophilic macromolecular chain-transfer agents capped with a trithiocarbonyl RAFT end group, which can be used as a steric stabilizer in the synthesis of well-defined amphiphilic block polymer nanoparticles via an additional RAFT polymerization step.

2.7.2 Preparation

Inform yourself about the hazard potential and the physical properties of the substances involved. Prepare a proper preliminary report. Prepare RAFT polymerization mechanism for the preliminary discussion. Try to find answers to the questions in the notes.

Chemicals:

4-acryloylmorpholine (NAM), azobisisobutyronitrile (AIBN), 4-cyano-4-[(dodecylsulfanylthiocarbonyl)sulfanyl]-pentanoic acid (CDPA), anisole, *N,N*-dimethylacetamide (DMAc), deuterated chloroform, diethyl ether.

2.7.3 Procedure

NAM (3 mmol, 100 equivalent), CDPA (0.03 mmol, 1 eq.), AIBN (0.006 mmol, 0.2 eq.) and anisole (0.125 mmol, internal standard) are dissolved in DMAc (3 mL) in a reaction tube.

Oxygen is removed by carefully bubbling nitrogen through the solution for 30 minutes. The reaction is started by placing it in a preheated oil bath at 70 °C and stirring for 2 hours. After stopping the polymerization by cooling in ice and air exposure, the mixture is analyzed via ^1H NMR spectroscopy in CDCl_3 to determine monomer conversion.

The product is purified via precipitation into diethyl ether to remove residual starting materials and is subsequently dried overnight in a vacuum oven. After drying, the molar mass of the polymer is determined by size-exclusion chromatography (SEC) in DMF and ^1H NMR spectroscopy. Additionally, the polymer solution is analyzed via UV-visible spectroscopy at 310 nm.

NMR samples are prepared by dissolving 1 drop of polymer solution or approx. 10 mg of purified polymer is dissolved in 0.5 mL CDCl_3 .

For SEC, the purified polymer is dissolved in DMF (final concentration of 1 mg/mL). After stirring 1 hour, the solution (1 to 2 mL) is passed through an SEC filter.

To prepare the UV-vis spectroscopy sample, 5 mg of purified polymer is dissolved in 2 mL of water.

2.7.4 Notes

- 1) Calculate the required amounts of starting materials.
- 2) The stock solutions of CDPA and AIBN are prepared in DMAc. Calculate the required amounts from the stock solutions and additional volume of DMAc needed to reach the same concentration. Show your calculations to the lab assistant.

3) Why is CDPA selected as a RAFT agent?

4) Why is AIBN amount lower than CDPA?

5) Why is anisole used? How is the conversion calculated?

6) The kinetics of RAFT polymerization are carried out by two students. For this, samples are taken at precise, predetermined time intervals (30 minutes, 1 hour, and 2 hours) to determine monomer conversion and molar mass of the polymer. What do you expect to observe?

7) Which molar mass average is obtained by ^1H NMR spectroscopy?


2.7.5 Questions that needs to be answered in the report

- Calculate the conversion and average molecular weight of obtained polymer.
- Show the results of the kinetics experiments in an $\ln(1/(1-\text{conversion}))$ vs time plot. What can you conclude from this plot?
- Why is 100 % conversion not aimed?
- Why are the molecular weights obtained by SEC and ^1H NMR spectroscopy different?
- Why is the polymer solution analyzed by UV-vis spectroscopy?

2.7.6 Checklist

^1H NMR spectroscopy	
Conversion (%)	
SEC	
UV-vis spectroscopy	

2.8 PNAM-*b*-PSty nanoparticles via polymerization-induced self-assembly (PISA) (only after 2.7 is completed)

2.8.1 Introduction and Aim

Block copolymers (BCPs) combine the properties of two different polymers and often lead to nanostructures through phase separation and self-assembly. Traditionally, the self-assembly step of amphiphilic BCPs is carried out in dilute solutions using co-solvent or solvent-switch methods. These are often carried out at low concentrations and require prior purification of the polymers. Polymerization-induced self-assembly (PISA) is a method of *in situ* nanoparticle synthesis in which formation of the second block simultaneously occurs, saving the time and eliminating the need for intermediate purification. It also enables much higher solids contents, up to around 50 wt%.

In PISA, a solvophilic precursor is extended with a monomer that forms a solvophobic block in a selective solvent, leading to self-assembly into various nanostructures such as spheres, rods/fibers, and vesicles. PISA is an advanced form of emulsion/dispersion polymerization and uses the reactive precursor block as a stabilizer, eliminating the need for traditional surfactants. The final morphology of the nanoparticles depends on parameters such as the block length ratio of solvophilic to solvophobic segments, solids content, and solvent selection, among others.

In this experiment, an aqueous RAFTPISA process will be used to form an amphiphilic block copolymer using the NAM-based, macroRAFT agent synthesized in the previous step as the hydrophilic precursor and styrene as the core-forming monomer.

2.8.2 Preparation

Inform yourself about the hazard potential and the physical properties of the substances involved. Prepare a proper preliminary report. Prepare the corresponding RAFT polymerization mechanism for the preliminary discussion. Try to find answers to the questions in the notes.

Chemicals:

PNAM-macroRAFT agent*, styrene (Sty), potassium peroxydisulfate (KPS), water, deuterated DMSO

* M_n of the polymer synthesized in the previous experiment needs to be calculated for the further calculations.

2.8.3 Procedure

For emulsion RAFTPISA of PNAM with styrene DP_{th} = 100, 10 wt% solids content

PNAM (0.016 mmol, 1 equivalent) and KPS (0.0032 mmol) are dissolved in required amount of water. Styrene (1.6 mmol, 100 eq.) is added and the mixture is stirred for 10 minutes to obtain a homogeneous emulsion.

After the solution is deoxygenated by nitrogen bubbling for 30 minutes, polymerization is carried out at 80 °C and 500 rpm for 24 hours and stopped by cooling down and air exposure.

The mixture is analyzed by ^1H NMR spectroscopy in deuterated DMSO to determine styrene conversion and the molar mass is determined via SEC in DMF (for the latter, the remaining of the PISA dispersion should be let dry overnight in the vacuum oven at 50 °C).

The polymer solution is analyzed via dynamic light scattering (DLS) at room temperature after dilution to 1 wt% concentration to determine the so-called Z-average of the nanoparticles.

The dried polymer is characterized by DSC.

The NMR sample is prepared by dissolving 1 drop of the reaction mixture in 0.5 mL DMSO-*d*₆.

The dried dispersion is dissolved in DMF (final concentration of 1 mg/mL). After stirring 1 hour, the solution is passed through an SEC filter.

For DSC measurement, 10 mg of dried dispersion should be weighed in a vial.

2.8.4 Notes

- 1) Calculate the required amounts of starting materials for DP_{Sty} 100, 200 and 400 at 10 wt%, and for DP_{Sty} 100 at 5 wt%, 10 wt% and 15 wt%. The water volume should be between 2.5 and 3 mL. Each student should carry out a different reaction condition and the results should be combined.
- 2) In this experiment, is it possible to carry out a random copolymerization instead of RAFT polymerization?
- 3) What is the difference between the traditional self-assembling methods of BCPs and PISA?
- 4) Why is the solvent selection important in a PISA process?
- 5) What is the Z-average and how does DLS roughly work?

2.8.5 Questions needs to be answered in the report

- What are the stages of emulsion RAFT PISA?
- Calculate the conversion of styrene and molar mass of the block copolymer.
- Prepare a table showing the reaction conditions and hydrodynamic diameters of the nanoparticles.
- What is the effect of hydrophobic block length on nanoparticle size? What is the effect of concentration on nanoparticle size? Discuss these effects in relation to your results.
- What can be T_g of block copolymer? What is the effect of high T_g on PISA process?

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